



Editorial Article

Journal of MAR Pulmonology (Volume 3 Issue 5)

## **Feasibility of Antigen Test Kits for Massive Screening and Detection of SARS-CoV-2 (COVID-19) Infection**

Attapon Cheepsattayakorn<sup>1,3\*</sup>, Ruangrong Cheepsattayakorn<sup>2</sup>, Porntep Siriwanarangsun<sup>3</sup>

1. 10th Zonal Tuberculosis and Chest Disease Center, Chiang Mai, Thailand
2. Department of Pathology, Faculty of Medicine, Chiang Mai University, Chiang Mai, Thailand
3. Faculty of Medicine, Western University, Pathumtani Province, Thailand

**Corresponding Author: Attapon Cheepsattayakorn**, 10th Zonal Tuberculosis and Chest Disease Center, 143 Sridornchai Road Changklan Muang Chiang Mai 50100 Thailand

**Copy Right:** © 2021 Attapon Cheepsattayakorn, this is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

**Received Date: September 13, 2021**

**Published date: October 01, 2021**

### **Introduction:**

N and S proteins of SARS-CoV-2 (COVID-19) are the main immunogens among the four structural proteins (E, M, N, and S) [1]. The N protein-IgG ELISA provides a sensitivity of 94.7 % that is significantly higher than that of the S protein-IgG ELISA [2]. Antibodies against N proteins are longer-lived and have greater volume, in comparison to E, M, and S proteins [2]. Approximately, 100 companies are manufacturing or developing rapid antigen tests (RATs), one of the four types (virus isolation with cell cultures, serological testing, RATs, and molecular techniques) [2] for SARS-CoV-2 (COVID-19) detection

Citation: Attapon Cheepsattayakorn. "Feasibility of Antigen Test Kits for Massive Screening and Detection of SARS-CoV-2 (COVID-19) Infection." MAR Pulmonology 3.5  
[www.medicalandresearch.com](http://www.medicalandresearch.com) (pg. 1)

[3]. Most RATs for SARS-CoV-2 (COVID-19) detection use a simple-to-use lateral flow test format that commonly used for influenza, malaria, and HIV testing as a sandwich immunodetection [4]. The testing results are interpreted by the operator within 10 to 30 minutes after collecting the respiratory specimen and applying it to the test strip [5]. In comparison to the nucleic acid amplification tests (NAATs), a decreasing sensitivity is found in the trade-off for simplicity of RATs for SARS-CoV-2 (COVID-19) operation [4]. As of September 11, 2020, only three companies submitted documents towards WHO's Emergency Use Listing (EUL) procedure, two tests have been approved by Japan's Pharmaceutical and Medical Devices Agency, and only four tests have received United States Food and Drug Administration (FDA) Emergency Use Authorization (EUA) [6, 7]. In comparison to NAATs in respiratory specimens (nasal or nasopharyngeal swabs), the specificity is consistently high (> 97 %), whereas the sensitivity ranges from 0 to 94 % [8-17]. When the cycle threshold (Ct) values are equal to or less than 25 or the SARS-CoV-2 (COVID-19) viral loads are more than 10<sup>6</sup> genomic virus copies/mL that frequently present in the pre-symptomatic period (1-3 days before the symptom onset) and the first 5-7 days of the acute COVID-19 illness phase [18-20]. Kweon et al conducted a study on evaluation of the diagnostic accuracy of the two newly developed, point-of-care, RATs, the ichroma COVID-19 AgTM and AFIAS COVID-19 Ag for detecting SARS-CoV-2 (COVID-19) infection by serially collecting 200 nasopharyngeal samples from 38 COVID-19-infected patients and 122 samples from negative control group [21]. The study revealed that both RATs demonstrated the sensitivity of 91.3 % to 100 % for samples with Ct < 25, whereas the specificity was 98.7 % to 98.9 % for AFIAS COVID-19 Ag and 100.0 % for ichroma COVID-19 AgTM [21]. The sensitivity of AFIAS COVID-19 Ag and ichroma COVID-19 AgTM for all targeted genes (E, N, and RdRP) was higher for samples collected before 7-days post-symptom onset than for those collected 8-14 post-symptom onset [21]. Stohr et al studied the sensitivity and specificity of the two self-testing kits (BD Veritor System-BD RAT, n = 1,604 and Roche SARS-CoV-2 RAT, n = 1,611) with lateral flow assay compared to the qtr.-PCR with Ct-value below the Ct-value cut-off demonstrated 78.0 % (95 % CI: 72.5-82.8) and 99.4 % (95 % CI: 99.0-99.6), respectively [22]. A test with the sensitivity of 80 % performed and implemented by at least 70 % of the population once a week was estimated to decrease the reproductive number of SARS-CoV-2 (COVID-19) from the basic reproductive number (R<sub>0</sub>) of 1.5 to the effective reproductive number (R<sub>e</sub>) below 1.0 [22].

The WHO has announced the general recommendations for the use of RATs for SARS-CoV-2 (COVID-19) detection as the following : 1) SARS-CoV-2 (COVID-19) RATs that meet the sensitivity of 80 % or higher and the specificity of 97 % or higher, compared to a NAAT reference assay where NAAT prolonged turnaround times preclude clinical utility or is unavailable; 2) Appropriate scenarios for the use of SARS-CoV-2 (COVID-19) RATs (responding to suspected COVID-19 outbreaks in remote settings, semi-closed communities, and institutions where NAAT is not immediately available, supporting outbreak investigations; testing of asymptomatic cases of contacts may be considered even if the RATs is not specifically authorized for this use, due to asymptomatic cases having been shown to have viral loads similar to symptomatic cases; where there is widespread community transmission; and monitoring

trends in COVID-19 incidence; 3) Initial introduction of RATs into clinical use; 4) in institutions where confirmed testing with NAAT is not feasible, and indications that RAT results may be incorrect should raise about validity suspicions; and 5) Use of RATs is not recommended in populations or settings with low expected prevalence of SARS-CoV-2 (COVID-19) (for examples; elective surgery, blood donation, screening at point of entry) [4]. The WHO also recommends the selection of RATs for procurement and implementation that includes 1) Quality of available data to validate the test; 2) Reported performance; 3) Manufacturing capacity and further evidence of quality; 4) Manufacturing quality and regulatory status; 5) Distribution and technical support; 6) Storage conditions, shelf-life, and shipping; 7) Availability, completeness, and clarity of instructions for use; 8) Cost of RATs; 9) Contents of RAT kit; and Specimen collection requirements [4]. Additionally, the WHO recommends the implementation considerations that include 1) Strictly following the supplier-recommended procedures; 2) Biosafety requirements for operators must be in place; 3) Each of these RATs has a specifically indicated method for specimen processing after collection; 4) Specimen collection is one of the most critical factors affecting performance of RATs; 5) Use of instrumented detection systems demands additional training requirements and sufficient infrastructure; and 6) Post-market surveillance with regulatory oversight [4]. Several variables may impact on RAT clinical performance that include 1) Pre-analytical influencers (sample type and way of sampling; collection device, transport media, and volume versus direct testing without dilution by transport media; time to test and storage/transport conditions, the time delay before processing); 2) Analytical influencers (Viral load of the specimen and viral load distribution in respective cohort (represented by SARS-CoV-2 (COVID-19) viral RNA copies/mL or Ct), analytical sensitivity and specificity of the PCR reference standard, PCR assay specifics, for examples, different targeted genes (E-/N-/RdRp-gene, etc., and across-laboratory differences (for examples, the definition of a positive specimen starting at Ct < 38 or < 40); and 3) Clinical parameters of the tested individual (days post-symptom onset of sampling, asymptomatic versus symptomatic status, the definition of symptoms “suspicion of SARS-CoV-2 (COVID-19) infection, and severity of symptom [23]. The sensitivity and specificity of the rapid SARS-CoV-2 (COVID-19) antigen (N protein) detection results are demonstrated in Table 1 [2].


<b>Country; Reference</b>	<b>Test; Manufacturer</b>	<b>Specimen Type</b>	<b>Method</b>	<b>Sensitivity; Specificity (%) (95 % CI)</b>
<b>Belgium; [24]</b>	<b>COVID-19 Ag Respi-Strip assay (Coris BioConcept, Gembloux, Belgium)</b>	<b>Nasopharyngeal secretions</b>	<b>Immunochromatographic test</b>	<b>30.2 (21.7- 39.9); 100</b>

Belgium; [25]	COVID-19 Ag Respi-Strip assay (Coris Bioconcept Gembloux, Belgium)	Nasopharyngeal secretions	Immunochromatographic test	30 (16.7-47.9); 100
Chile; [8]	SARS-CoV-2 antigen test (Bioeasy Biotechnology Co., Shenzhen, China)	Nasopharyngeal, oropharyngeal swabs	Fluorescence immunochromatographic test	93.9 (86.5-97.4); 100 (92.1-100)
France; [26]	COVID-19 Ag Respi-Strip assay (Coris BioConcept, Gembloux, Belgium)	Nasopharyngeal secretions	Immunochromatographic test	50.0 (39.5-60.5); 100 (91.8-100)
Germany; [27]	SARS-CoV-2 rapid antigen test (Roche, Switzerland)	Nasopharyngeal swabs	Lateral flow assay	70.7 (59.0-81.0); 96 (89.0-99.0)
USA; [28]	Sofia 2 SARS antigen FIA (Quidel Corporation)	Nasopharyngeal swabs	Lateral flow immunofluorescent sandwich assay	80 (68-88); Not available

**Table 1:** Demonstrating the sensitivity and specificity of the rapid antigen (N protein) detection tests for SARS-CoV-2 (COVID-19)

Recently, the United States Centers for Disease Control and Prevention established an interim guidance for RATs for detecting SARS-Cov-2 (COVID-19) infection in a community setting [29, Figure 1]. If the individual has a low risk of SARS-CoV-2 (COVID-19) infection, a positive RAT result from an asymptomatic individual in a community setting may need confirmatory testing (for examples; a low risk of SARS-CoV-2 (COVID-19) infection would be an individual who is fully COVID-19 vaccinated, or has had a SARS-CoV-2 (COVID-19) infection in the last 3 months, or has had no known exposure to a COVID-19-infected individual within the last 14 days) [29]. If a asymptomatic individual has a high risk

of SARS-CoV-2 (COVID-19) infection, a negative RAT result may need confirmatory testing (for examples; a high risk of SARS-CoV-2 (COVID-19) infection would be an individual who is not fully COVID-19 vaccinated, or has not had a SARS-CoV-2 (COVID-19) infection in the last 3 months, or has had close contact or suspected exposure to SARS-CoV-2 (COVID-19) within the last 14 days) and should consider performing serial RATs every 3-7 days for 14 days [29]. A suspected COVID-19-symptomatic individual with a positive RAT result can be interpreted that the individual is infected with SARS-CoV-2 (COVID-19) [29]. For fully COVID-19 vaccinated individual with a positive RAT, the healthcare providers should inform the public health authorities and ideally, a separate specimen would be collected and sent to a laboratory for public-health-purpose viral sequencing [29]. A positive RAT result in a symptomatic individual with a low risk of SARS-CoV-2 (COVID-19) infection (for examples; for examples; a low risk of SARS-CoV-2 (COVID-19) infection would be an individual who is fully COVID-19 vaccinated, or has had a SARS-CoV-2 (COVID-19) infection in the last 3 months, or has had no known exposure to a COVID-19-infected individual within the last 14 days) may need confirmatory testing [29]. A symptomatic individual with a negative RAT result should be performed serial RATs every 3-7 days for 14 days, in addition to confirmation with a laboratory-based NAAT [29]. In a symptomatic individual with a negative RAT result, the negative RAT result may not need confirmatory SARS-CoV-2 (COVID-19) testing [29]. If a symptomatic individual with a negative RAT result had not any known exposure to SARS-CoV-2 (COVID-19), the individual do not need to quarantine as well as for asymptomatic individual with a negative RAT result who had been fully COVID-19 vaccinated, or had a SARS-CoV-2 (COVID-19) in the last 3 months, or have has close contact or suspected exposure to a COVID-19-infected individual within the last 14 days [29].



**Figure 1:** Demonstrating RATs algorithm for community setting

(**Source:** United States Centers for Disease Control and Prevention. Interim guidance for antigen testing for SARS-CoV-2. Updated on May 13, 2021. Available at: <https://www.cdc.gov/lab/antigen-tests-guidelines> (accessed on July 31, 2021))

Currently, which SARS-CoV-2 (COVID-19) antigens are most appropriate for serological testing and whether the direct detection of viral antigens in the clinical specimens is a more sensitive-rapid-accurate-immunological diagnostic method for SARS-CoV-2 (COVID-19) remain unknown. S1 subunit protein of SARS-CoV-2 (COVID-19) is the S protein as an antigen for the serological diagnosis of SARS-CoV-2 (COVID-19) infection, whereas the S2 subunit protein plays a significant role in the cross-reactivity when the whole S protein is utilized as an antigen. Combinations of both S and N proteins have revealed improvement of the laboratory results through multiantigen protein arrays in comparison to each individual protein antigen. RATs become useful diagnostic tool for the SARS-CoV-2 (COVID-19)-early detection.

In conclusion, they should be used in conjunction with the molecular methods and further urgent studies should be focused on strategies to improve the accuracy and sensitivity and post-implementation evaluation of the diagnostic accuracy.

## **Bibliography**

1. Meyer B, Drosten C, Müller MA. Serological assays for emerging coronaviruses: challenges and pitfalls. *Virus Res* 2014; 194: 175-183. DOI: <https://doi.org/10.1016/j.viruses.2014.03.018>
2. Li D, Li J. Immunologic testing for SARS-CoV-2 infection from the antigen perspective. *Journal of Clinical Microbiology* 2021; 59 (5): e02160-20. DOI: <https://doi.org/10.1128/JCM.02160-20>
3. Foundation for Innovative New Diagnostics. SARS-CoV-2 diagnostics pipeline 2020. Available at <https://www.finddx.org/covid-19/pipeline/> (accessed on July 31, 2021).
4. World Health Organization. Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays: interim guidance. Published on September 11, 2020. Available at: <https://www.who.int/publications/items> (accessed on July 31, 2021).
5. Peto T, UK COVID-19 Lateral Flow Oversight Team. COVID-19: rapid antigen detection for SARS-CoV-2 by lateral flow assay: a national systematic evaluation of sensitivity and specificity for mass-testing. *E-Clinical Medicine* 2021; 36: 100924. DOI: <https://doi.org/10.1016/j.eclinm.2021.100924> Available at: <https://www.journals.elsevier.com/eclinicalmedicine> (accessed on July 31, 2021).
6. U.S. Food and Drug Administration. In vitro diagnostics EUAs 2020. Available at: <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/vitro-diagnostics-euas> (accessed on July 31, 2021).

7. Agency PaMD. PMDA's effort to combat COVID-19 2020. Available at: <https://www.pmda.go.jp/english/about-pmda/0002.html> (accessed on July 31, 2021).
8. Porte L, Legarraga P, Vollrath V, Aguilera X, Munita JM, Araos R, et al. "Evaluation of novel antigen-based rapid detection test for the diagnosis of SARS-CoV-2 in respiratory samples". *Int J Infect Dis* 2020; 99: 328-333. S1201-9712(20)30405-7.
9. Diao B, Wen K, Chen J, Liu Y, Yuan Z, Han C, et al. "Diagnosis of acute respiratory syndrome coronavirus 2 infection by detection of nucleocapsid protein". *medRxiv*. 2020: 2020.03.07.20032524 (accessed on July 31, 2021).
10. Lambert-Niclot S, Cuffel A, Le Pape S, Vauloup-Fellous C, Morand-Joubert L, Roque-Afonso AM, et al. "Evaluation of a rapid diagnostic assay for detection of SARS-CoV-2 antigen in nasopharyngeal swabs". *J Clin Microbiol* 2020; 58 (8).
11. Mertens P, De Vos N, Martiny D, Jassoy C, Mirazimi A, Cuypers L, et al. "Development and potential usefulness of the COVID-19 Ag Respi-Strip Diagnostic Assay in a pandemic context". *Front Med (Lausanne)* 2020; 7: 225.
12. Blairon L, Mokrane S, Wilmet A, Dessilly G, Kabamba-Mukadi B, Beukinga I, et al. "Large-scale, molecular and serological SARS-CoV-2 screening of healthcare workers in a 4-site public hospital in Belgium after COVID-19 outbreak". *J Infect* 2020. S0163-4453(20)30514-4.
13. Mak GC, Cheng PK, Lau SS, Wong KK, Lau CS, Lam ET, et al. "Evaluation of rapid antigen test for detection of SARS-CoV-2 virus. *Journal of Clinical Virology: the official publication of the Pan American Society for Clinical Virology* 2020"; 129: 104500.
14. Nagura-Ikeda M, Imai K, Tabata S, Miyoshi K, Murahara N, Mizuno T, et al. "Clinical evaluation of self-collected saliva by Rt-qPCR, direct RT-qPCR, RT-LAMP, and a rapid antigen test to diagnose COVID-19. *J Clin Microbiol* 2020". JCM.01438-20.
15. Omi K, Takeda Y, Mori M. "SARS-CoV-2 qRT-PCR Ct value distribution in Japan and possible utility of rapid antigen testing kit". *medRxiv*. 2020. 2020.06.16.20131243
16. Scohy A, Anantharajah A, Bodéus M, Kabamba-Mukadi B, Verroken A, Rodriguez-Villalobos H. "Low performance of rapid antigen detection test as frontline testing for COVID-19 diagnosis". *J Clin Virol* 2020; 129 : 104455.
17. Weitzel T, Legarraga P, Iruretagoyena M, Pizarro G, Vollrath V, Araos R, et al. "Head-to-head comparison of four antigen-based rapid detection tests for the diagnosis of SARS-CoV-2 in respiratory samples". *bioRxiv*. 2020. 2020.05.27.119255

18. Weiss A, Jellingso M, Sommer MOA. "Spatial and temporal dynamics of SARS-CoV-2 in COVID-19 patients: a systematic review and meta-analysis". *EBioMedicine* 2020; 58.
19. Arons MM, Hatfield KM, Reddy SC, Kimball A, James A, Jacobs JR, et al. "Pre-symptomatic SARS-CoV-2 infections and transmission in a skilled nursing facility". *N Engl J Med* 2020; 382 (22) : 2081-2090.
20. Dinnes J, Deeks JJ, Adriano A, Berhane S, Davenport C, Dittrich S, et al. "Rapid, point-of-care antigen and molecular-based tests for diagnosis of SARS-CoV-2 infection". *Cochrane Database of Systematic Reviews* 2020 (8).
21. Kweon OJ, Lim YK, Kim HR, Choi Y, Kim M-C, Choi S-H, et al. "Evaluation of rapid SARS-CoV-2 antigen tests, AFIAS COVID-19 Ag and ichroma COVID-19 Ag, with serial nasopharyngeal specimens from COVID-19 patients". *PLoS ONE* 2021; 16 (4) : e0249972.
22. Stohr JJJM, Zwart VF, Goderski G, Meijer A, Nagel-Imming CRS, Kluytmans-van den Bergh MFQ, et al. "Self-testing for the detection of SARS-CoV-2 infection with rapid antigen tests. medRxiv preprint. 2021. Posted on February 23, 2021. DOI: <https://doi.org/10.1101/2021.02.21.21252153> Available at: <https://www.medrxiv.org/content> (accessed on July 31, 2021).
23. Hayer J, Kasapic D, Zemmrch C. Real-world clinical performance of commercial SARS-CoV-2 rapid antigen tests in suspected COVID-19: a systematic meta-analysis of available data as of November 20, 2020. *International Journal of Infectious Diseases* 2021; 108: 592-602. DOI: <https://doi.org/10.1016/j.ijid.2021.05.029> Available at: <https://www.ijidonline.com/article/fulltext>, [www.elsevier.com/locate/ijid](http://www.elsevier.com/locate/ijid) (accessed on July 31, 2021).
24. Scohy A, Anantharajah A, Bodé M, Kabamba-Mukadi B, Verroken A, Rodriguez-Villalobos H. "Low performance of rapid antigen detection test as frontline testing for COVID-19 diagnosis". *J Clin Virol* 2020; 129 : 104455. DOI : <https://doi.org/10.1016/j.jcv.2020.104455>
25. Blairon L, Wilmet A, Beukinga I, Tré-Hardy M. "Implementation of rapid SARS-CoV-2 antigenic testing in a laboratory without access to molecular methods : experiences of a general hospital". *J Clin Virol* 2020; 129 : 104472. DOI : <https://doi.org/10.1016/j.jcv.2020.104472>
26. Lambert-Niclot S, Cuffel A, Le Pape S, Vaulou-Fellous C, Morand-Joubert L, Roque-Afonso AM, et al. "Evaluation of a rapid diagnostic assay for detection of SARS-CoV-2 antigen in nasopharyngeal swab". *J Clin Microbiol*; 58 : e00977-20. DOI : <https://doi.org/10.1128/JCM.00977-20>
27. Krüttgen A, Cornelissen CG, Dreher M, Hornef MW, Imohl M, Kleines M. "Comparison of the SARS-CoV-2 rapid antigen detection test to the Real Star SARS-CoV-2 RT-PCR kit". *J Virol Methods* 2020; 288 : 114024. DOI : <https://doi.org/10.1016/j.jviromet.2020.114024>

28. United States Food and Drug Administration. Individual EUAs for antigen diagnostic tests for SARS-CoV-2. Available at : <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/vitro-diagnostics-euas#individual-antigen> (accessed on July 31, 2021).

29. United States Centers for Disease Control and Prevention. Interim guidance for antigen testing for SARS-CoV-2. Updated on May 13, 2021. Available at : <https://www.cdc.gov/lab/antigen-tests-guidelines> (accessed on July 31, 2021).